



## User Guide

### HFS1000 Non-Contact Infrared Thermometer

HFS1000 is included as a medical device on the Australian Register of Therapeutic Goods (ARTG No. 279076).

HFS1000 is intended for use as a measuring device for determining the patient's body temperature. **ALWAYS READ THE LABEL AND USER GUIDE, USE ONLY AS DIRECTED.**

Keep this user guide for future reference. Always retain your proof of purchase in case of warranty service and register your product on line at AUSTRALIA: [www.oricom.com.au](http://www.oricom.com.au)





## **Need Help?**

If you need assistance setting up or using your Oricom product now or in the future, call Oricom Support.

Australia      1300 889 785 or 02 4574 8888  
[www.oricom.com.au](http://www.oricom.com.au)  
Mon-Fri 8am – 6pm AEST

New Zealand    0800 67 42 66  
[www.oricom.co.nz](http://www.oricom.co.nz)  
Mon-Fri 10am – 8pm NZST

**Table of Contents**

---

HOW TO TAKE A TEMPERATURE	10
PRODUCT DESCRIPTION	14
BEFORE FIRST USE	16
SWITCHING ON AND SETTING THE THERMOMETER	16
MEASURING BODY TEMPERATURE	17
DISPLAYING STORED MEASUREMENTS	20

## 4 **HFS1000** NON-CONTACT INFRARED THERMOMETER

MEASURING OBJECT TEMPERATURE	21
MEASURING ROOM TEMPERATURE	23
ERROR MESSAGES	24
BATTERY REPLACEMENT	26
MAINTENANCE AND STORAGE	27
PRODUCT SPECIFICATIONS	28
CUSTOMER SUPPORT	31

# **HFS1000 NON-CONTACT INFRARED THERMOMETER** **5**

IMPORTANT	31
HOW TO MAKE A CLAIM UNDER WARRANTY IN AUSTRALIA	32
IMPORTANT INFORMATION	33
WARRANTY INFORMATION (AUSTRALIA)	33

Thank you for purchasing the Oricom HFS1000 Non-contact Infrared Thermometer.

This product has been developed using the latest in infrared technology. This technology allows the user to measure the temperature without physical contact. To use this device correctly please follow the directions in this user guide. The instructions provide important details on how to achieve accurate measurements.

## Warnings




Please read this user guide carefully before initial use.

- Use this product only for its intended purpose as described in this manual and follow all instructions.
- The measurement of temperature using this product does not represent a doctors medical examination.
- Consult with your medical professional if any abnormality is found.
- Do not take a temperature measurement over a scar or open tissue.
- Do not expose the product to temperatures below 20 °C (68 °F) or above 50 °C (122 °F) or above humidity of 95% ( $\geq 95\%Rh$ ).
- Do not drop the thermometer as it may cause damage.
- Do not take the temperature with this product near heat sources such as in direct sunlight, near a fireplace or stove.

- Do not use this product if it is not working properly, contact Oricom support.
- There are no user serviceable parts the only service required is battery replacement and general cleaning please follow instructions in this user guide.
- If the product is not going to be used for an extended period of time, please remove the batteries and store the thermometer in a cool dry place.
- This product is NOT waterproof. Direct contact with water or other liquids should be avoided.
- Do not leave the thermometer with infants or children at any time.
- The protective glass over the lens is fragile, handle with care. Do not touch the glass of the infrared lens with bare fingers.
- Do not modify or repair the thermometer.
- Many people may not have an elevated temperature even if they are ill. In the following cases, please consult your medical professional before using this product.
  - > Under 90-day-old infants.
  - > People on steroids, antibiotics or antipyretics.
  - > People with compromised immune systems.

## Cautions

- The reliability of the measurement cannot be guaranteed if the temperature is measured from parts of body that are not recommended.
- Similar to other thermometers, measurement error may occur within a range of  $\pm 0.2\text{ }^{\circ}\text{C}$  ( $\pm 0.4\text{ }^{\circ}\text{F}$ ).
- Measurement should not be taken until completely relaxed at normal room temperature.
- Please note that measurement may be incorrect after exercising, bathing or getting wet in the rain.
- Temperature between right and left temple could be different.
- Measurement may be incorrect when low battery () signal is indicated. Replace the battery and re-measure.
- When measuring the temperature of a liquid (milk, water), place the probe as close as possible to the liquid but ensure that it is not immersed into the liquid as this may damage your thermometer.
- Check your product regularly for damage and exhausted batteries.

### NOTE:

ASTM laboratory accuracy requirements in the display range of 22 to 40  $^{\circ}\text{C}$  (71.6  $^{\circ}\text{F}$  to 104.0  $^{\circ}\text{F}$ ) for IR thermometers is  $\pm 0.3\text{ }^{\circ}\text{C}$  ( $\pm 0.5\text{ }^{\circ}\text{F}$ ), whereas for mercury-in-glass and electronic thermometers, the requirement per ASTM Standards E 667-86 and E 1112-86 is  $\pm 0.1\text{ }^{\circ}\text{C}$  ( $\pm 0.2\text{ }^{\circ}\text{F}$ ).



## Simple Operation Principle

The body temperature is reflected on the temporal artery.

- The temporal artery is the artery located closest to the skin and is spread across the forehead. If the blood with body temperature reflection passes the temporal artery on forehead, corresponding amount of infrared rays will be generated. From the infrared rays, the body temperature can be measured.
- This product is a contactless type thermometer and does not make direct contact with forehead.
- The temperature measurement can be taken 2 to 3 cm away from the body or object.
- This will decrease discomfort or possible infections caused by touch.

# 10 HFS1000 NON-CONTACT INFRARED THERMOMETER

## How to take a temperature

### IMPORTANT

#### MEASUREMENT POINT

- Recommended measurement is the temple area.



Temple area

Or



Behind the ear lobe



Incorrect measurement



2~3 cm

- ▶ 2~3 cm away from forehead.



Wipe away sweat



Push back hair

- ▶ If sweating on the forehead is present we recommend changing the measurement point to behind to earlobe. Wipe away sweat on the forehead and push back any hair from the forehead before use.



Behind the ear lobe

- ▶ If you have applied an ice or cool pack, you need to wait for at least 30 minutes before taking a temperature measurement.
- ▶ The measurement should then be taken behind the ear lobe.



- ▶ Do not measure temperature near heat or cooling sources.  
Do not use outdoors.



WAIT 30 Minutes

- ▶ If there is a significant change of surrounding temperature, please, allow the thermometer to adjust to the new temperature for at least 30 minutes before use.


## 12 HFS1000 NON-CONTACT INFRARED THERMOMETER

- Body Temperature Measurement Ranges

Body temperature means internal temperature of body. Normal temperature even in a healthy person can vary within a small range. The temperature of a person measured in the evening is normally 0.5 °C higher than early morning and is also affected by the environment, time and activity. Refer to the table below to find out the normal body temperature ranges.

### Normal ranges by site:

• Axillary:	94.5 – 99.1 °F	34.7 – 37.3 °C	• Rectal:	97.9 – 100.4 °F	36.6 – 38.0 °C
• Oral:	95.9 – 99.5 °F	35.5 – 37.5 °C	• HFS1000:	96.8 – 100.2 °F	36.0 – 37.9 °C



- Using the Thermometer
  - **Measuring Body Temperature**  
Select the HFS1000 to body temperature measuring mode, and check the body temperature on the temple area from 2 to 3 cm distance.
  - **Measuring Bathing Water Temperature**  
Select the thermometer to daily wide temperature measuring mode, and check the water temperature by holding the thermometer 2 to 3 cm from the water. Be careful not to get water on the thermometer.
  - **Measuring the Nursing Bottle**  
Select the thermometer to daily wide temperature measuring mode, and check the temperature by holding the thermometer 2 to 3 cm from the nursing bottle.


**CAUTION:**

ALWAYS CHECK BATH WATER OR NURSING BOTTLE TEMPERATURE MANUALLY IN ADDITION TO USING THIS THERMOMETER.

# 14 HFS1000 NON-CONTACT INFRARED THERMOMETER

## Product Description

Unit



- 1 LCD display
- 2 Sound (on/off) button
- 3 Mode button
- 4 Memory button
- 5 Start button (Power & Measurement)
- 6 Strap position
- 7 Probe

### Buttons Functions



Switching device on. Starting temperature measurement.



Setting the sound (On/Off).




Selecting the measurement method and unit.



Displaying stored measurements (max 32).  
Entering the measurement method and unit.

Display



- 1 Body temperature mode
- 2 Object temperature mode
- 3 Room temperature mode
- 4 Low Battery symbol
- 5 Temperature display
- 6 Measurement unit (Celsius/Fahrenheit)
- 7 Memory
- 8 Sound on/off

# 16 HFS1000 NON-CONTACT INFRARED THERMOMETER





## Before first use

---

- Before first use, insert the batteries supplied.  
> Battery Replacement [page. 26]

## Switching on and setting the thermometer

---

- Press the Start button  or Mode button  briefly.
  - > The device switches on. After two seconds, the display shows that the device is ready for use.
  - > The device always starts in body temperature mode. 
  - > After 1 minute without use, device automatically turns off.  
How to force the device to power off: Press and hold  for 5 seconds.







## Setting the measurement unit (°C/°F)

- Press  3 times to select Centigrade and 4 times to select Fahrenheit.  
If correct measurement unit (°C or °F) is selected, press the MEM button  to save it.






## Setting the measurement sound

- Press the Sound button  briefly.
  - > Each time the button is pressed, the sound switches on or off.
  - > If the display shows  after pressing the button, the device emits sound.
  - > If  is not displayed after pressing the button, the device does not emit sounds.
  - > The device always starts in "Sound on" mode. 


## Measuring body temperature

---


- Press the Start button  or Mode button  briefly.
  - > The device switches on. After two seconds, the display shows that the device is ready for use.
  - > The device always starts in body temperature mode. 
- Hold the thermometer 2 to 3 cm from the measuring point.
  - > The ideal measuring point is the temple area.



# 18 HFS1000 NON-CONTACT INFRARED THERMOMETER

- Press the Start button  briefly.
  - > Keep the distance until you hear the measurement completion beep.
  - > If the thermometer is in silent mode, keep the distance for approximately 2 seconds.

## Normal



**Site of measurement recommended is the temple area**

## Sweating



**behind the ear lobe**

### NOTE:

If sweat is present, we recommend changing the measurement point from the temple to behind the ear lobe as this area is less affected by sweating.

- Temperature readings may vary according to location. Therefore we recommend to measure the temperature at more than 2 points. Then please recognize the highest temperature as the right temperature.
- Check the temperature on the LCD screen. The temperature reading will be automatically stored.
- The LED backlight will illuminate based on the readings.
  - > LED will automatically turned off after 2 to 3 seconds.

## LED 3 Colours

Green: between 22 °C ~ 37.3 °C

Yellow: between 37.4 °C ~ 37.9 °C







Red: higher than 38 °C  
lower than 22 °C




## 20 HFS1000 NON-CONTACT INFRARED THERMOMETER

### Displaying stored measurements

---

- Press the Start button  or Mode button  briefly to turn on the thermometer.
  - > The device switches on. After two seconds, the display shows that the device is ready for use.
- The device automatically stores the last 32 measurements.
- Press the MEM button  briefly.
  - > The most recent measurement is displayed.
- Each additional press of the MEM button  calls up the next most recent measurement. The number on the display indicates how many measurements ago this reading was taken.
- Briefly press the Start button  or Mode button  to exit this feature.
  - > The device switches on. After two seconds, the display shows that the device is ready for use.
- The device only stores measurements when in body temperature mode.



## Measuring object temperature









- Press the Start button  or Mode button  briefly.
  - > The device switches on. After two seconds, the display shows that the device is ready for use.
  - > The device always starts in body temperature mode. 
- Press the Mode button  briefly.
  - > The device switches to the object temperature measuring mode. 
  - > If the object temperature measuring mode  is selected, press the MEM button  to save it.
- Check the measurement ready status in [Figure 1].
- Hold the thermometer 2 to 3 cm from the desired measuring point.
- Press the Start button  briefly.
  - > Keep the distance until you hear the measurement completion beep.
  - > If the thermometer is in silent mode, keep the distance for approximately 2 seconds.



Figure 1



## 22 HFS1000 NON-CONTACT INFRARED THERMOMETER

- Read the temperature from the display. [Figure 2]
  - > Measurements in object temperature mode are not stored.

### Object temperature






- **LED 3 Colours**

- > Green: Measured temperature is between 0 °C ~ 100 °C (32 °F ~ 212 °F)
- > Red: Measured temperature is higher than 100 °C (212 °F)
- > Red: Measured temperature is lower than 0 °C (32 °F)



Figure 2

## Measuring Room temperature

- If you want to use the thermometer to measure the room temperature, you need to change the mode to room temperature. 
- Press the Mode button  briefly.
  - > The device switches to the room temperature measuring mode. 
  - > If the object temperature measuring mode  is selected, press the MEM button  to save it.
- Check the measurement ready status in [Figure 1].
- Read the temperature from the display [Figure 2].
  - > Measurements in room temperature mode are not stored.

### Room temperature

- **LED 3 Colours**

- > Green: Measured temperature is between 16 °C ~ 40 °C (60.8 °F ~ 104 °F)
- > Red: Measured temperature is higher than 40 °C (104 °F)
- > Red: Measured temperature is lower than 16 °C (60.8 °F)




Figure 1







Figure 2

## 24 HFS1000 NON-CONTACT INFRARED THERMOMETER

### Error messages

Display	Error description	Trouble-shooting
	<p>Ambient temperature is not within the allowed operating temperature.</p> <p>Operating temperature range: 16 °C ~ 40 °C (60.8 °F ~ 104 °F)</p> <p>“H 40” means: Ambient temperature is 43 °C (109.4 °F) which is higher than allowed range.</p> <p>“L 16” means: Ambient temperature is 12 °C (53.6 °F) which is lower than allowed range.</p>	<p>Allow the device to remain for approx 30 minutes in a normal room temperature.</p>




Display	Error description	Trouble-shooting
	The temperature of the measuring point is above the measurement range.	Measurement range for body temperature: 22 °C UP 42.5 °C (71.6 °F UP 108.5 °F) Measurement range for object temperature: 0 °C UP 100 °C (32 °F UP 212 °F)
	The temperature of the measuring point is below the measurement range.	Measurement range for body temperature: 22 °C UP 42.5 °C (71.6 °F UP 108.5 °F) Measurement range for object temperature: 0 °C UP 100 °C (32 °F UP 212 °F)
  <div>Figure 1</div> <div>Figure 2</div>	Figure 1. The batteries are weak.	Change batteries. Figure 2. Display appears briefly after changing the batteries.

## 26 HFS1000 NON-CONTACT INFRARED THERMOMETER

### Battery Replacement

- The device requires two AA Alkaline batteries.
  - > Open the battery compartment.
  - > For this purpose, slide down the battery compartment lid on the back.
- Take the used batteries out of the battery compartment. Insert new batteries.
  - > Observe correct polarity.
- Close battery compartment.
  - > Used batteries should not be disposed of in normal household waste.
  - > Dispose of batteries according to local recycling regulations.



Important, if equipment is not likely to be used for some time.

- > Storage and transfer conditions, temperature: -20 °C to +50 °C
- > Remove the batteries if the device is not to be used for over three months, the batteries may leak and cause a malfunction.

## Maintenance and Storage

---

- The sensor part of the thermometer is very sensitive and must be kept clean and undamaged.
- Do not put any foreign materials inside the probe.
- To clean the sensor, using a cotton swab with alcohol (ethanol), carefully wipe the surface clean.
- Keep the product away from direct sunlight, dust and any other type of contaminating substance.
- The thermometer should be stored in normal room temperature.
- Please place the product in normal room temperature for 30 minutes prior to use if it has been stored outside normal room temperature.



Do not subject the thermometer to strong shocks, such as dropping the unit.  
Do not submerge the device or any of the components in water or any other solutions.  
Store the thermometer and the components in a clean, safe location.  
Do not modify this thermometer without authorization of the manufacturer. Do not disassemble or attempt to repair the unit or components.

# 28 HFS1000 NON-CONTACT INFRARED THERMOMETER

## Product Specifications

Item	Details of description
Measurement method	Infrared (IR) measurement, Non-Contact measurement
Basic functions	Switches for Multi-measurements (Body, Object and Room temperature) Switches for Celsius/Fahrenheit, Memory function, Sound controls (On,Off)
Operating condition	16 °C to 40 °C (60.8 °F to 104 °F), Humidity: 15% to 93%, non-condensing, Atmospheric Pressure (hPa) 700~1060
Transport and Storage condition	-20 °C to 50 °C (-4 °F to 122 °F), Humidity: 15% to 93%,non-condensing Atmospheric Pressure (hPa) 700~1060
Measuring distance	2 to 3 cm from the measuring point is the temple area
Measuring range	<ul style="list-style-type: none"><li>• Body temperature: 22 °C to 42.5 °C (71.6 °F to 108.5 °F), Accuracy: <math>\pm 0.2</math> °C/<math>\pm 0.4</math> °F (36 °C~39 °C/196.8 °F~102.2 °F), Other Accuracy: <math>\pm 0.3</math> °C/<math>\pm 0.5</math> °F</li><li>• Object temperature: 0 °C to 100 °C (32 °F to 212 °F), Accuracy: <math>\pm 3</math> °C/<math>\pm 5.4</math> °F</li><li>• Room temperature: 16 °C to 40 °C (60.8 °F to 104 °F), Accuracy: <math>\pm 3</math> °C/<math>\pm 5.4</math> °F</li></ul>

# HFS1000 NON-CONTACT INFRARED THERMOMETER 29

Item	Details of description
Measuring time	Less than 2 seconds, Continuous measurement
Resolution	0.1 °C/0.1 °F
Mode selection	Body temperature, Object temperature, Room temperature
Memory function	Automatic memory function, Confirm the result by using memory button Can store up to 32 results (1 to 32)
Sound	Button sound, Measuring sound, Result sound, The sound can be turned on and off by using Sound button
Screen type	Liquid crystal display (mono LCD)
Display	Body temperature mode, Object temperature mode, Room temperature mode, Memory location, Low Battery, °C/°F degree, Hi / Lo signs for the temperature outside the measurable ranges

## 30 HFS1000 NON-CONTACT INFRARED THERMOMETER

Item	Details of description
Button	Start, Sound, Mode, Memory
Backlight	Turns On for 3 seconds after measurement. The backlight turns off automatically
Power-saving	Turns Off automatically after one minute
Components	HFS1000 1 set, Manual, Battery (AA type) 2ea, case
Power	=== AA Alkaline Batteries (x2)
Measurable numbers	More than 10,000 times

## Customer Support

---

If you have any problems setting up or using this product you will find useful tips and information in this user guide as well as “Frequently Asked Questions” on our website [www.oricom.com.au](http://www.oricom.com.au).

If you have further questions about using the product after reviewing the resources above or would like to purchase replacement parts or accessories please call our Customer Support Team. Our dedicated local support team are more likely to be able to help you than the retailer where you made your purchase.

## Important

---

Please retain your purchase receipt and attach to the back page of this user guide as you will need to produce this if warranty service is required. Take a few moments to register your product online: [www.oricom.com.au](http://www.oricom.com.au)

### How to make a claim under Warranty in Australia

---

Oricom has a simple warranty process for you to follow:

- Please call or email our Customer Support Team, contact details follow.
- A Customer Support Team member will verify after troubleshooting with you if your product qualifies under warranty. If so, they will give you a Product Return Authorisation number.
- We will then email or fax a Return Authorisation form and a Repair Notice (if necessary), together with instructions on how to return the goods for warranty service.

Please note that if a Customer Support Team member advises that your product does not qualify for return, this warranty does not apply to your product.

Products that are authorised to be returned to Oricom in Australia must include all of the following:

- A completed Return Authorisation form
- A copy of your Proof of Purchase (please keep your original copy)
- The faulty product, including all accessories.

Send the approved returns to:

Oricom International Pty Ltd  
Locked Bag 658  
South Windsor NSW 2756 Australia



Please note that this warranty excludes expenses incurred by you in returning any faulty product to us. You must arrange and pay any expenses incurred (including postage, delivery, freight, transportation or insurance of the product) to return the faulty product to us, however, we will arrange delivery of the repaired or replaced faulty product to you.

## **Important Information**

---

### **Repair Notice**

Please be aware that the repair of your goods may result in the loss of any user-generated data (such as stored telephone numbers, text messages and contact information). Please ensure that you have made a copy of any data saved on your goods before sending for repair.

Please also be aware that goods presented for repair may be replaced by refurbished goods or parts of the same type rather than being repaired.

### **Warranty Information (Australia)**

---

This Warranty is provided by Oricom International Pty Ltd ABN 46 086 116 369, Unit 1, 4 Sovereign Place, South Windsor NSW 2756, herein after referred to as "Oricom".

## 34 **HFS1000 NON-CONTACT INFRARED THERMOMETER**

Oricom makes no other warranties or conditions, express or implied, including as to acceptable quality and fitness for a particular purpose, except as stated in this Warranty.

Any implied warranties that may be imposed by law are limited in duration to the Warranty Period.

Oricom warrants that the product is free from defects in materials or workmanship during the Warranty Period. This Warranty does not extend to any product from which the serial number has been removed or was purchased outside of Australia.

This warranty in no way affects your statutory warranty rights under the Competition and Consumer Act 2010 or any other similar legislation.

The Warranty Period will be 1 year from the date of purchase of the product evidenced by your dated sales receipt. You are required to provide proof of purchase as a condition of receiving warranty services.

You are entitled to a replacement product or repair of the product according to the terms and conditions of this document if your product is found to be faulty within the Warranty Period. This Warranty extends to the original purchaser only and is not transferable.

Products distributed by Oricom are manufactured using new materials or new and used materials equivalent to new in performance and reliability. Spare parts may be new or equivalent to new. Spare parts are warranted to be free from defects in material or workmanship for thirty (30) days or for the remainder of the Warranty Period of the Oricom branded product in which they are installed, whichever is longer. During the Warranty Period, Oricom will where possible repair and if not replace the faulty product or part thereof. All component parts removed under this

Warranty become the property of Oricom. In the unlikely event that your Oricom product has a recurring failure, Oricom may, subject to the Competition and Consumer Act 2010, at its discretion, elect to provide you with a replacement product of its choosing that is at least equivalent to your product in performance.

Oricom does not warrant that the operation of the product will be uninterrupted or error free.

Oricom is not responsible for damage that occurs as a result of your failure to follow the instructions that came with the product. These terms and conditions together with any specific terms and conditions contained in the user guide to the product purchased constitute the complete and exclusive agreement between you and Oricom regarding the product.

No change to the conditions of this Warranty is valid unless it is made in writing and signed by an authorised representative of Oricom.

Oricom will not be in breach of a warranty expressly set out in this User Guide or under the Competition and Consumer Act 2010 and excludes any liability for damages or any other remedy arising under any other legislation or the common law if the damage occurs as a result of:

1. failure by you to adhere to the warnings and follow the instructions set out in this user guide for the proper installation and use of the product;
2. negligence on your part or misuse by you of the product;
3. an uncontrollable external cause which results in the product not functioning including but not limited to power failure, lightning or over voltage; and

## 36 HFS1000 NON-CONTACT INFRARED THERMOMETER

4. modification to the product or services carried out on the production by anyone other than Oricom or Oricom's authorised service provider.

Oricom will not be liable for any damages caused by the product or the failure of the product to perform, including any lost profits or savings or special, incidental or consequential damages. Oricom is not liable for any claim made by a third party or made by you on behalf of a third party. This limitation of liability applies whether damages are sought, or a claim made, under this Warranty or as a tort claim (including negligence and strict product liability), a contract claim or any other claim. However, this limitation of liability will not apply to claims for personal injury. Nothing in this Warranty excludes, restricts or modifies any condition, warranty, right or remedy which pursuant to the Competition and Consumer Act 2010 applies to this Warranty and which may not be so excluded, restricted or modified. For warranties that cannot be excluded, restricted or modified, Oricom limits the remedies available to those specified in the relevant legislation.

Oricom products come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.



## **Contact details for Oricom support and warranty claims in Australia**

Oricom International Pty Ltd  
Locked Bag 658  
South Windsor, NSW 2756  
Australia

Email: [support@oricom.com.au](mailto:support@oricom.com.au)  
Phone: 1300 889 785 or (02) 4574 8888  
(Monday to Friday 8am to 6pm AEST)  
Web: [www.oricom.com.au](http://www.oricom.com.au)  
Fax: (02) 4574 8898

## **Contact details for Oricom support and warranty claims in New Zealand**

Email: [support@oricom.co.nz](mailto:support@oricom.co.nz)  
Phone: 0800 674 266  
(Monday to Friday 10am to 8pm NZST)  
Web: [www.oricom.co.nz](http://www.oricom.co.nz)



Ref: 21092016