


oricom

IET800

INFRARED

EAR THERMOMETER



COLOUR
CODED
READING
BASED
ON AGE

The IET800 is included as a medical device on the Australian Register of Therapeutic Goods (ARTG 525670).



PRE-HEATED
PROBE



NIGHT
MODE



360° FLOATING SELF
ADAPTION PROBE

oricom

IET800

INFRARED

EAR THERMOMETER



KEY FEATURES

0+ SUITABLE FOR ALL AGES
Gentle enough for newborns, designed to grow with your baby.

COLOUR CODED READING BASED ON AGE

Age groups	Normal	Moderate fever	High fever
	1 beep	2 beeps	3 beeps
0 - 3 months	35.8°C ≤ T ≤ 37.4°C	/	T > 37.4°C
3 - 36 months	35.4°C ≤ T ≤ 37.6°C	37.6°C < T ≤ 38.5°C	T > 38.5°C
Above 36 months	35.4°C ≤ T ≤ 37.7°C	37.7°C < T ≤ 39.4°C	T > 39.4°C

PRE-HEATED PROBE
Warms to body temperature before entering the ear.

360° FLOATING SELF-ADAPTION PROBE
Self-adapts to any ear canal shape.

DISPOSABLE PROBE COVER WITH EJECTOR
One-touch ejector removes used probe covers without contact.

ERGONOMIC AND COMFORTABLE DESIGN
Comfortable to hold and easy to manoeuvre.

FAST ACCURATE TEMPERATURE DETECTION (APPROX 3-5 SECS)*
Accurate temperature readings in seconds.

LARGE BACKLIT LCD DISPLAY

STORES UP TO 32 RECORDED TEMPERATURES
Monitor fever progress and changes at a glance.

AUTOMATIC POWER OFF AFTER APPROX 60 SECONDS OF NO ACTIVITY

INCLUDES PROTECTIVE CASE

REQUIRES 2 AA ALKALINE BATTERIES (INCLUDED)

INCLUDES PACK OF 20 PROBE COVERS

12 MONTH WARRANTY
Backed by a 12 month warranty for added peace of mind.

PACK INCLUDES



1 X EAR THERMOMETER

1 X PROTECTIVE CASE

20 X DISPOSABLE PROBE COVERS

2 X AA ALKALINE BATTERIES

1 X USER MANUAL

REFILLS



40 X PROBE COVERS (IETPC)

This product is intended for use as a measuring device for determining the patient's body temperature. **ALWAYS READ THE LABEL AND USER GUIDE, USE ONLY AS DIRECTED.**
CAUTION: This product contains a lithium-ion battery, see user guide for details.

This product is manufactured in People's Republic of China under a quality system by: Guangzhou Berrcom Medical Device Co., Ltd, No. 38, Huanzhen Xi Road, Dagang Town, Nansha, 511470, Guangzhou, People's Republic of China.

IMPORTED & DISTRIBUTED IN AUSTRALIA AND NEW ZEALAND by the sponsor:-
Oricom International Pty Ltd, Unit 1, 4 Sovereign Place, South Windsor 2756 AUSTRALIA. Phone: 02 4574 8888

oricom

Oricom® is the registered trade mark of Oricom Holdings Pty Ltd. All other trade marks are trade marks of their respective owners. Due to continued product development, the product may differ slightly to that illustrated in this brochure.

oricom.com.au